



## Centri-Sifter® High Temperature Sterilize-in-Place Unit

**KASON CORPORATION**  
67-71 East Willow St.  
Millburn, NJ 07041-1416  
USA  
Tel: 973-467-8140  
Fax: 973-258-9533  
E-mail: [info@kason.com](mailto:info@kason.com)

**KASON CORPORATION,  
EUROPE**  
Block 4, Units 12 & 13  
Park Hall Business Village  
Park Hall Road  
Longton, Stoke-on-Trent  
Staffordshire, ST3 5XA  
UNITED KINGDOM  
Tel: 1782-597540  
Fax: 1782-597549  
E-mail:  
[sales@kasoneurope.co.uk](mailto:sales@kasoneurope.co.uk)

**SEPARATOR ENGINEERING LTD.**  
6025 Boul Thimens  
Saint Laurent  
Quebec  
CANADA  
H4S 1V8  
Toll Free (Canada only):  
877-694-4441  
Tel: 514-667-6777  
Fax: 514-745-2074  
E-mail:  
[info@separatorengineering.com](mailto:info@separatorengineering.com)



[www.kason.com](http://www.kason.com)

# CENTRI-SIFTER® High Temperature Sterilize-in-Place Unit

This ASME coded sanitary CENTRI-SIFTER separator de-lumps and scalps fine pharmaceutical powders. It is designed to provide 300°F (149°C) temperatures and 60 psig of internal pressure required for in-place sterilization. Using special materials and finishes enables this unit to meet the specifications for pharmaceutical, biotechnical vaccine applications.

## Options:

- Anti-blinding devices
- Bag Dump Station
- Clean-in-place (CIP)
- Twin centrifugal configuration
- Design, construction and finish to FDA, 3-A, BISSC, EEC and other sanitary standards
- Belt drive and Direct drive configurations

